


















 <p>BATH BENCH W/BACK ALUM LEG 300LB CAPACITY ADJ SEAT HT: 15.5-19.5" 10320-4, \$49.95</p>	 <p>BATH BENCH ALUM LEG 300LB CAPACITY ADJ SEAT HT: 15.5-19.5" 10323-4, \$32.50</p>	<p>Guardian</p>  <p>BATH BENCH ADJ WITH BACK, PADDED ARMS 350LB CAP, ADJ HT 16-21" MDS89745RAH, \$65.50</p>	<p>Guardian</p>  <p>BATH BENCH ADJ, PAD ARMS 350LB CAP, ADJ HT 16-21" MDS89740RWA, \$49.95</p>
 <p>SHOWER STOOL ROUND SEAT ALUMINUM LEG ADJUST 1" INCREMENTS 300LB CAP 10324-4, \$29.95</p>	 <p>SHOWER STOOL SWIVEL SEAT WITH LOWER TRAY ADJ HT 18-20", 300LB CAP 1042, \$79.95</p>	 <p>BATH BENCH HD WITH BACK 500LB CAPACITY ADJ HEIGHT: 14-19" 10321-4, \$79.95</p>	 <p>BATH BENCH HD 500LB CAPACITY ADJ HEIGHT: 14-19" 10325-4, \$69.95</p>
 <p>SHOWER CHAIR 360° SWIVEL SEAT ARMS+BACK; BUILT-IN ANTIMICROBIAL AGENT 4-POS LOCK, 350LB CAP RTL12A001-GR, \$169.95</p>	 <p>TRANSFER BENCH MAXI-DRN 350LB CAP HT ADJ: 18-1/4-22-1/4" 7927KD-4, \$85.50</p>	 <p>TRANSFER BENCH ECON OPEN SEAT, 400LB CAP HT ADJ: 17.75-21.75" 12011KDC-1, \$129.95</p>	 <p>TRANSFER BENCH PADDED GREY 300LB CAP HT ADJ: 17.5-22.5" 7955KD-1, \$149.95</p>
 <p>TRANSFER BENCH SLIDING SEAT FOLD/LEGS ADJ 21-25", 300LB CAP RTL12075, \$185.95</p>	 <p>TRANSFER BENCH WITH SLIDING SEAT AND FOLD UP LEGS drive</p>	 <p>TRANSFER BENCH HD 500LB CAPACITY HT ADJ 18.5-23.5" 10328-4, \$149.95</p>	

 <p>ELEV TOILET SEAT ELONG W/REMOV ARMS LOCK 5" HEIGHT 300LB CAP 12013, \$57.50</p>	 <p>ELEVATED TOILET SEAT ECONOMY ADDS 5.5" HEIGHT B302, \$24.50</p>	 <p>ELEV TOILET SEAT LOCKING SEAT HT: 4.5", 300LB CAP B305, \$49.95</p>	 <p>ELEVATED TOILET SEAT LOCK WIDER OPENING ADDS 4.25" HEIGHT, 500LB CAP B31300, \$59.95</p>
 <p>ELEVATED TOILET LOCKING WITH ADJ WIDTH ARMS ADDS 5" HEIGHT FITS MOST TOILETS, 300LB CAP B30300, \$75.50</p>	 <p>ELEV TOILET SEAT W/ARMS ELONG, 250LB CAP USES EXISTING SEAT&LID, ADDS 3.5" HEIGHT B31600, \$75.50</p>	 <p>ELEV TOILET SEAT 3.5"H, INSTALLS UNDER SEAT/LID, 300LB CAP B307 (FOR ROUND BOWLS) B306 (ELONGATED) \$39.95</p>	 <p>ELEV TOILET SEAT+ARMS INSTALLS UNDER SEAT/LID, 300LB CAPACITY 1043 (ROUND), 1044 (ELONG), \$75.95</p>
 <p>ELEVATED TOILET SEAT WITH ARMS, LOCK, 600LB CAPACITY 55CBAS0028R, \$79.95</p>	 <p>TOILET SEAT CUSHION VELCRO FASTEN TO SEAT 1040, 2" THICK, \$25.95 1041, 4" THICK, \$29.95</p>	 <p>TOILET SAFETY FRAME 18-25"WIDE, ADJ HT: 28-32" 10352-4, \$39.95</p>	 <p>TOILET SAFETY FRAME COMPACT 16-18" WIDE 300LB CAPACITY 10351-4, \$39.95</p>
 <p>COMMODE SAFETY RAIL PADDED ARM MAGAZINE HOLDER LVA1023, \$69.95</p>	 <p>TUB RAIL CLAMP-ON HT ADJ: 14.5-17" 300LB CAP RTL12036-ADJ \$49.95</p>	 <p>GRAB BAR FLIP UP SHURGRIP ALL STAINLESS STEEL WITH HARDWARE 97-2230G, \$239.95</p>	 <p>THE TRAVELER 4" SUCTION MOUNT HANDLE 1168, \$15.95</p>

	40033A, 12", \$19.50		2012A, 12", \$17.95
	40433A, 16", \$19.95		2016A, 16", \$18.50
	40133A, 18", \$21.00		2018A, 18", \$18.95
	40233A, 24", \$22.95		2019A, 24", \$21.95
	40333A, 32", \$22.95		2021A, 32", \$22.95
<p>HEAVY DUTY STEEL GRAB BARS WITH WHITE ENAMEL FINISH 250LB CAP, 1.5" FROM WALL</p>		<p>CHROME PLATED HEAVY DUTY STEEL GRAB BARS WITH KNURLED FINISH, LIFETIME WARRANTY 250LB CAP, 1.5" FROM WALL</p>	

			
<p>GRAB BAR SUCTION MOUNT ULTRA WITH INDICATOR - CHANGES RED TO GREEN TO TELL WHEN PROPERLY FASTENED 12" LONG, B20000, \$18.95 16" LONG, B19600, \$19.95</p>		<p>WALL MOUNT SUCTION GRAB BARS 12": 1067, \$22.95 16": 1167, \$26.95</p>	<p>SHOWER HOLDER EASY- RELEASE TAB 1039, \$14.95</p>

			<table border="1"> <thead> <tr> <th colspan="2">SPECIFICATIONS</th> </tr> </thead> <tbody> <tr> <td>Width Between Arms</td> <td>17 1/2"</td> </tr> <tr> <td>Overall Width</td> <td>24 1/2"</td> </tr> <tr> <td>Seat Width</td> <td>14"</td> </tr> <tr> <td>Seat Depth</td> <td>15"</td> </tr> <tr> <td>Seat Height Range</td> <td>15.5" - 23.5"</td> </tr> <tr> <td>Weight Capacity</td> <td>300 LB</td> </tr> </tbody> </table>	SPECIFICATIONS		Width Between Arms	17 1/2"	Overall Width	24 1/2"	Seat Width	14"	Seat Depth	15"	Seat Height Range	15.5" - 23.5"	Weight Capacity	300 LB
SPECIFICATIONS																	
Width Between Arms	17 1/2"																
Overall Width	24 1/2"																
Seat Width	14"																
Seat Depth	15"																
Seat Height Range	15.5" - 23.5"																
Weight Capacity	300 LB																
<p>HAND HELD SHOWER SPRAY 84" WITH ON-OFF BUTTON B21500, \$29.95</p>	<p>HAND HELD SHOWER 84" WITH ON/OFF SWITCH 5 SETTINGS 9311-R, \$29.95</p>	<p>DURABLE GREY-COATED STEEL FRAME COMPLETE WITH PAIL & COVER AND SPLASH SHIELD 3-IN-1 COMMODE 7103A-4, \$69.95</p>															

	<ul style="list-style-type: none"> • Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame • Welded steel construction provides strength and stability • Increased seat depth provides added stability and comfort • Durable, white-powder coated steel frame is attractive and easy to clean • Contoured waterfall armrests are designed to enhance user comfort and provide a good gripping surface 								
	<table border="1"> <tr> <td>Width Between Arms</td> <td>18 1/2"</td> </tr> <tr> <td>Overall Width</td> <td>24 1/2"</td> </tr> <tr> <td>Overall Depth</td> <td>20"</td> </tr> <tr> <td>Seat Width</td> <td>14"</td> </tr> </table>	Width Between Arms	18 1/2"	Overall Width	24 1/2"	Overall Depth	20"	Seat Width	14"
	Width Between Arms	18 1/2"							
	Overall Width	24 1/2"							
	Overall Depth	20"							
Seat Width	14"								
<table border="1"> <tr> <td>Seat Depth</td> <td>16 1/2"</td> </tr> <tr> <td>Seat Height Range</td> <td>17" - 27"</td> </tr> <tr> <td>Weight Capacity</td> <td>300 LB</td> </tr> </table>	Seat Depth	16 1/2"	Seat Height Range	17" - 27"	Weight Capacity	300 LB			
Seat Depth	16 1/2"								
Seat Height Range	17" - 27"								
Weight Capacity	300 LB								
<p>3-IN-1 COMMODE DLXE 7103, \$79.95</p>									



COMMODE FOLDING STEEL
350LB CAP ADJ HT: 15.75-21.5"
10371, \$69.95



COMMODE DROP-ARM 3-IN1 STEEL
300LB CAPACITY ADJ HT: 20-24"
11125KD-2, \$129.95



Folding frame design assembles in seconds without tools
Seat Height: 18 1/4" to 23 1/2"
Width between the arms: 24"
Overall width: 28"

COMMODE BARIATRIC FOLDING
650LB WT CAPACITY
G1201LBX1, \$169.95



The Dynarex Bariatric Drop Arm Bedside Commode is a versatile, three-in-one commode that can be configured as a commode, raised toilet seat or toilet safety rail over your current toilet. The tool-free design allows for easy set up with height adjustments in 1-inch increments with drop arms allowing for safe lateral patient transfers. This durable, steel welded construction includes a plastic snap on seat with lid, 10-quart drop-in bucket with handle, splash guard and removable tool-free back.

- Constructed from heavy-duty, powder-coated steel
- Easy-release drop-arm mechanism allows for lateral patient transfers
- Commode bucket slides in from the front of the commode for added convenience
- Includes commode bucket with handle, cover, and splash shield
- Tool-free setup and height adjustment
- Weight Capacity 1000 lb.
- Limited Lifetime Warranty

10375. \$225.00

Cleanis CareBag® Commode Liner with Super Absorbent Pad



The single-use CareBag® Commode Liner with Super Absorbent Pad is designed to safely collect, transport, and dispose of body fluids, thus avoiding unpleasant hand cleaning of the pail. The super absorbent pad turns liquid waste (urine, diarrhea, blood, serous fluid) into a gel within seconds. When used to line the commode, both the user and caregiver are protected from the spread of pathogenic organisms by direct or indirect contact, thus minimizing risk of cross infection. 20.5" L x 15-7/10" W, Polyethylene, disposable. Eliminates odors when completely closed
7831738, \$21.00 -BX/20



Commode Pail Liners 19x11.5"
Holds up to 2 quarts of liquid,
tie handles
P70900, \$11.95 -BX/7