

Description	Add'l Info	Item #	Price
-------------	------------	--------	-------

**PROSPHYG 768 POCKET ANEROID SPHYG with MATCHING CASE**



**The 768 series features:**

- Contemporary black enamel 300mmHg manometer provides day-in, day-out dependability exceeds industry standards
- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system, prevents mis-cuffing
- Chrome-plated brass air release valve
- Extra large, color-matched, dual pocket carry case, measures 10" x 6" x 3"
- Not made with natural rubber latex
- Assembled, inspected, and packaged in the U.S.A. from foreign components

- 768-11ARB..... ROYAL BLUE
- 768-11ADG..... HUNTER GREEN
- 768-11AN..... NAVY
- 768-11AV..... PURPLE
- 768-11ABK..... BLACK
- 768-11AR..... RED
- 768-11ABD..... BURGUNDY
- 768-11AOR..... ORANGE

**\$55.95**

ANEROID SPHYG BLACK NYLON CUFF LF	LG ADULT 34-50CM	ADC 720-12XBK	79.95
ANEROID SPHYG ADIMALS NYLON CUFF LF	ADULT 23-40CM	ADC 720-11AAD	69.95
ANEROID SPHYG B/C AWARE NYLON CUFF LF	ADULT 23-40CM	ADC 720-11ABCA	69.95
ANEROID SPHYG BLACK NYLON CUFF LF	ADULT 23-40CM	ADC 760-11ABK	45.95
ANEROID SPHYG BLACK NYLON CUFF LF	CHILD 13-19.5CM	ADC 760-9CBK	45.95

**A&D AUTOMATIC BLOOD PRESSURE MONITORS**



**Product Features**

- Clinically Validated for Accuracy
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- One Button Operation
- Store 60 Readings each for up to 4 users
- Average reading feature provides trending
- Time & date stamp for every reading
- AccuFit Plus Cuff fits most arms 8.6-16.5"
- Multi-user , upper arm

**PREMIUM WRIST BLOOD PRESSURE MONITOR**

Advanced, multi-user blood pressure monitor that provides trusted accuracy and reliable measurements every time **UB543, \$79.95**

**UA767FAC (WITH AC ADAPTER).....\$89.95**  
**UA767F (WITHOUT AC ADAPTER).....\$79.95**

DIGITAL BP MONITOR AUTO XL CUFF 16.5-23.6"	WITH MEMORY,AC ADAPTER	UA789AC	189.95
DIGITAL BP MONITOR AUTO SMALL ADULT CUFF	WITH MEMORY,AC ADAPTER	UA767PSAC	89.95



OTOSCOPE,POCKET XENON 3X MAG	INCLUDES:3-SPECS,2-AAA BAT	PRESTIGE 3470	65.00
PENLIGHT WITH PUPIL GUAGE DISP	EACH,6/PK	ADC 351	2.95
PENLIGHT REUSE METALITE WHT PUPIL	2-AAA BAT INCLUDED	ADC 352WP	6.95

**3470 XENON POCKET OTOSCOPE**

**PULSE OXIMETERS**



- 2 Year Warranty
- Free Case Included!
- Easy to Read Large LCD Display
- FDA Approved for Accuracy
- Efficient - Over 30 hours (1 AAA Battery Required)



**DB12, \$89.95**



Fingertip measurement of SpO2% (Oxygen Saturation) and pulse rate. Measures quickly and accurately. LED display screen can display information in two directions. Low power consumption with auto shut-off feature. Up to 30 hours of continuous use. Lightweight and compact design.

- 459-BLK BLACK
- 459-HKW HOT PINK
- 459-PUR PURPLE
- 459-REW RED/WHITE

**\$39.95**

Description	Add'l Info	Item #	Price
<b>STETHOSCOPES</b>			
 <p><b>LITTMANN CARDIOLOGY - Diagnose and make decisions with confidence</b></p> <ul style="list-style-type: none"> <li>• Hear subtle changes in patient status</li> <li>• Useful in critical care and challenging environments such as the ED, ICU, Cardiac ICU, step-down unit and other dynamic locations</li> <li>• Ergonomic and high profile construction; 40% larger chestpiece and 60% deeper bell than the 3M™ Littmann® Classic III™ Stethoscope</li> <li>• Designed for use with adult and pediatric patients</li> <li>• Tunable, dual-sided stainless steel chestpiece with open or closed bell</li> <li>• Non-chill bell sleeve for greater patient comfort</li> <li>• Soft-sealing eartips provide optimal comfort and excellent sound occlusion</li> <li>• Made in the USA without natural rubber latex</li> </ul>			
STETHOSCOPE,LITTMANN CARDIOLOGY 4	NAVY TUBE 27"	3M 6154	279.95
STETHOSCOPE,LITTMANN MASTER CARDIOLOGY	BLACK TUBE 27"	3M 2160	265.35
 <p>Littmann's latest addition to their line. <b>The Littmann Classic III</b> features patented tunable technology that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece. High acoustic sensitivity for performing general physical assessments. The bell stays clear of dirt and debris by keeping it covered with the small diaphragm.</p>			
STETHOSCOPE LITTMANN CLASSIC 3	BLACK TUBE 27"	3M 5620	126.50
STETHOSCOPE LITTMANN CLASSIC 3	NAVY TUBE 27"	3M 5622	126.50
STETHOSCOPE LITTMANN CLASSIC 3	CARIBBEAN BLUE TUBE 27"	3M 5623	126.50
STETHOSCOPE LITTMANN CLASSIC 3	BURGUNDY TUBE 27"	3M 5627	126.50
STETHOSCOPE LITTMANN CLASSIC 3	CIEL BLUE TUBE 27"	3M 5630	126.50
STETHOSCOPE LITTMANN CLASSIC 3	LAVENDER TUBE 27"	3M 5632	126.50
STETHOSCOPE LITTMANN LITEWEIGHT BLK	BLACK TUBE 28"	3M 2450	60.95
STETHOSCOPE LITTMANN LITEWEIGHT BRG	BURGUNDY TUBE 28"	3M 2451	60.95
STETHOSCOPE LITTMANN LITEWEIGHT CEIL	CEIL BLUE TUBE 28"	3M 2454	60.95
STETHOSCOPE LITTMANN LITEWEIGHT CRB	CARRIBEAN BLUE TUBE 28"	3M 2452	60.95
STETHOSCOPE LITTMANN LITEWEIGHT PINK	PINK TUBE 28"	3M 2456	60.95
STETHOSCOPE CLINICAL PLUS LTW BLK	31",4.9oz SS STEM	S123-BLK	39.95
STETHOSCOPE CLINICAL PLUS LTW PINK	31",4.9oz SS STEM	S123-HPK	39.95
STETHOSCOPE CLINICAL PLUS LTW NAVY	31",4.9oz SS STEM	S123-NAV	39.95
STETHOSCOPE CLINICAL PLUS LTW PURPLE	31",4.9oz SS STEM	S123-PUR	39.95

**The Adtemp™ 418N features:**

- 8-second measurement using proprietary predictive technology
- Range 93.2°F to 111.0°F ±.2°F or 34.0°C to 43.9°C ±.1°C depending upon scale selection
- Large LCD display with color-coded backlighting (green for readings in the normal range and red for readings above 99.5°F/37.5°C)
- Replaceable 3V battery (CR2032) provides up to 2,000 measurements
- Measures 5.31"L x 1.14"W x 0.59"H, weighs 0.77 oz (22 g) with battery
- Integral scabbard-type carry case, Includes five probe sheaths
- Meets the accuracy requirements for intermittent thermometers under ASTM E1112

**ADC 418N, \$16.95**



**NON-CONTACT THERMOMETERS**




- Non-Contact, Fast , Accurate Readings
  - Measures Body, Surface and Ambient Temperatures
  - Backlit Digital LCD Display Screen
  - High Temperature Alarm
  - Measures in Both Fahrenheit & Celsius
  - Features 32 data memory
- Includes two AAA batteries, 1 year warranty

**DT-29, \$59.95**



- 2 modes, 1 push: Body & object mode
- 1-second accurate measurement time
- Auto turn-off function for energy saving
- 32 readings memory recall 1 push switch
- Fahrenheit °F or Celsius °C

**SM017011CF, \$39.95 -Each**

Description	Add'l Info	Item #	Price
 <p>The Adtemp™ 429 Non-Contact Thermometer features an infrared sensor that reads temperature from the middle of the forehead with no skin contact required for a safe and comfortable experience that minimizes the risk of cross-contamination.</p> <p><b>The Adtemp™ 429 Infrared Thermometer features:</b></p> <ul style="list-style-type: none"> <li>• Infrared sensor simplifies temperature measurement—ideal for screening</li> <li>• Measures in about 1 second without patient contact eliminating the concern of cross contamination or costly probe covers</li> <li>• Measurement feedback system signals proper operating distance for accuracy</li> <li>• Measures at a distance of 4-6cm (1.5"-2.5") from patient's forehead</li> <li>• Dual mode - body and object</li> <li>• Large backlit LCD display is easy to read</li> <li>• Wide Temp Range</li> <li>• 30 reading memory</li> <li>• Fahrenheit-Celsius switchable</li> <li>• Auto-off function after 30 seconds conserves battery life</li> <li>• Audible beep when measurement is complete, Night mode mutes volume</li> <li>• Compact ergonomic design 5" x 2.06" x 1.75" (129.9mm x 52.4mm x 44.5mm), weighs approximately 3oz (87g) with batteries</li> <li>• Optional cradle available for wall mounting</li> </ul> <p>Uses 2 AAA batteries (included)</p> <p><b>ADC 429, \$ 85.95</b>                      <b>2 year warranty</b></p>			