



FOLDING WALKER WITH WHEELS, 350LB CAPACITY
10210-4 ADULT: ADJ HT 32-39"
10211-4 YOUTH: ADJ HT 25-32"
\$75.95



HD FOLDING WALKER WITH WHEELS 600LB CAPACITY
604070W, \$125.00



HEMI WALKER FOLDING
ADJ HT 30-35" 300LB CAP
6015A, \$59.95



TRAVEL WALKER RED/BLACK
ADJUST HEIGHT, 300LB CAP
10167, \$79.95



WALKER CLEVERLITE 5"
CASTERS 300LB CAP
ADJ SEAT HT: 20.5-24.5"
10230, \$159.95



KNEE WALKER STEERABLE HAND
BRAKES 300LB CAP
10216, \$225.00



WALKER TRAY WITH CUP
HOLDER 23"Wx1.5"Hx17"D
10124, \$27.95



WALKER BASKET
WITH INSERT
438B, \$19.95



WALKER CADDY
15.5"Wx5"Hx12.5"D
RTL10131, \$29.95



WALKER PLATFORM
ATTACHMENT
HT ADJ: 37-53"
6132B, \$69.95



WALKER BALLS GRAY,
RED, BLUE, PURPLE
\$8.99 -PR



WALKER SKI
GLIDES
3010, \$19.95 -PR



WALKER HAND GRIPS
SUP1064, \$16.95 -PR

- ✓ Contoured foam padding
- ✓ Breathable, moisture-wicking materials
- ✓ Adjustable for a secure fit
- ✓ Easy to clean
- ✓ Latex-free

Made in USA



The Lumex Set n' Go Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one. Designed with multiple height settings on the frame, the Set n' Go allows the seat height to be adjusted from a hemi 18" to a standard 22" height in 1" increments. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated. Adjustable handle height with knobs to secure in place, easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock). 6" casters, 300lb. weight capacity. Available in burgundy, blue, silver and pink
RJ4700, \$129.95



The Lumex Walkabout Four-Wheel Contour Deluxe rollator offers a comfortable, high end and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 16.5 lbs and supports up to a 300 lb. weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport. Ergonomic hand grips, easy to operate- locking and "quick stop" loop brakes (pull up to stop; push down to lock), padded seat for added comfort. 8" casters, contour back. Available in champagne, black, blue, burgundy, pink, lavender
RJ4805, \$159.95



- Fold-Up and Removable Back Support, New seamless padded seat
 - Zippered pouch under seat for added privacy and security of personal items
 - Removable, hinged, padded backrest can be folded up or down
 - 6" Black non-marring casters with soft grip tires are ideal for indoor and outdoor
 - Easy-to-use deluxe loop locks
 - Brakes with serrated edges provide firm hold
- R726 colors: Black, Blue, Red, Silver \$129.95**

UPRIGHT ROLLATOR WALKER

- Safely walk upright
- Easy foldup design
- Adjustable handles
- Comfortable seating
- Removable shoulder bag

Armrest height: 41-46"
Overall width: 25"
Seat width: 17.5"
Wheels: 8" diameter
Weight capacity: 300lbs
Lightweight aluminum frame
Total weight: 20lbs
MOB1033, \$259.95





- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage handle is located on seat
- Plastic clip holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- Handsome and removable zippered storage bag. Unique bag attachment keeps bag securely in place when rollator is open or folded
- Lightweight, aluminum frame
- Extra-tall frame models available

NITRO ROLLATOR

RTL1026, Available in red, white and black. 300lb Weight Capacity, **\$325.00**



Empower Rollator

- Padded backrest flips forward and backward enabling the user to sit comfortably on either side.
- Seat flips up to allow user to stand or walk inside the rollator
- Folds completely flat for easy storage and transport
- Features easy-to-identify yellow adjustment points, folding cup holder, memory foam seat and neoprene storage bag

MDS86845B, \$225.00



Guardian Deluxe Rollators

- Deluxe padded seat & curved padded backrest provide enhanced support and comfort
- Large 8" (20 cm) wheels combined with lightweight aluminum frame allows for smooth operation on most terrain.
- 300-lb weight capacity
- 17.5" (44 cm) distance between handles; seat height is 20" (51 cm); min./max. height 31-36.5" (79-93 cm). Weighs 20.5 lbs. (9 kg) Available in Blue or Rose

G07887, \$179.95



Simplicity 2 Rollator

- Rollator features adjustable handle height and backrest height for a custom fit
- Simply push down and lock brakes to use as a seat
- Smooth-rolling wheels make mobility a breeze
- Comes with bag and backrest

MDS86890EUB, BLUE
MDS86890EUS, SILVER

\$199.95

Imperial 4-Wheel Rollator HD

- Lightweight design makes rollator easy to use, fold and store for convenient transport
- Padded seat and backrest provide additional comfort
- Contoured backrest designed to fit the curvature of a user's back
- Ergonomic hand grips and hand brakes allow user to operate with ease
- Complete with pouch for storage of personal items
- Large 8" wheels are ideal for indoor and outdoor use
- 500LB weight capacity
- Available in various red, blue and black
- Newly added reflectors increase visibility

RJ4405, \$295.00



Basic Aluminum Rollators

- Easy-to-use loop-style brakes; push down to lock, pull up to release
- Padded backrest and seat for comfortable resting
- Height-adjustable arms for a custom fit
- Folding design for compact storage
- Available in Burgundy, Blue, Black, Purple, Green

MDS86850E, \$119.95



3-Wheel Rollator

Comes standard with basket, tray and pouch. Basket conveniently folds for easy storage. Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob. 7.5" casters for indoor or outdoor use. Handles easily adjust in height with self-threading knob. Special loop lock made of internal aluminum casting operates easily and ensures safety. Easy one-hand folding Available in red or blue

10289, \$129.95

Use As A Rollator
Or Wheelchair



ROLLATOR HYBRID – Use as a rollator or wheelchair

- Versatile hybrid rollator wheelchair allows individuals to walk, ride and rest safely on any terrain, including grass, gravel, carpeting and more
- Features removable footrests and adjustable handles for a comfortable fit. Adjustable height handles 33.25" – 38.25"
- Constructed with a durable, corrosion-resistant aluminum, the combination rollator transport chair safely supports up to 300 lbs.
- Folding frame easily fits in the trunk of a car and in compact storage areas

MOB1018, \$279.95

WOODEN CANES



7/8"X36" MOHAGANY

7/8"X36" BLACK

1"X36" NATURAL

1"X36" WALNUT

1"X36" WALNUT STEPPED/SCORCH

1"X36" ROSEWOOD

1"X42" EXTRA LONG WALNUT, CHESTNUT

LADIES 3/4" BLACK, ROSE

ALUMINUM ADJUSTABLE CANES



HURRYCANE FREEDOM EDITION FOLDING CANE, \$49.95

- 8 convenient height adjustments
- 360 degree pivoting head
- Freestanding design
- Anodized aluminum shaft easily folds for storage
- Wrist strap and travel bag included
- 350lb cap, adjustable height 30.5-37.5"
- Available in black, blue, red and purple



9021908 9051408 9051658 9051908
9052115 9052121 9052133 9052142

9021908, Black Alum Finger grip handle
Adjustable 34-42", \$25.50

9051408, Black Ergon Grip, Adjustable 30-39" \$29.95

Aluminum Canes Hypalon Grip & Strap adjustable 30-39", \$29.95

9051658, Black/Granite
9052115, Floral
9052121, Dogwood
9052127, Black/Marble
9052133, Paisley
9052142, Wild Flowers
9052188, Mauve

OFFSET HANDLE CANES



American Military & Pride Cane, \$19.95

LUMEX 5970CA Camouflage

LUMEX 5970FL American Flag

(Height adjustable 30-39" 250 lb weight capacity)



5940A 6220A 6221A 6327

ALUMINUM CANES ADJUSTABLE HEIGHT

5940A, 31-39" Offset Grip \$18.50

6220A, 30-39" Round Grip \$18.50

6221A, 30-39" Round Soft Grip \$18.50

6327, 31-39" Ortho-Ease grip bronze finish, \$25.20

QUAD CANES



9202508 Quad Cane
Small Base Adj 29-38"
Black, \$48.67

9202608 Quad Cane
Large Base Adj 30-39"
Black, \$35.50



6121A Quad Cane
Small Base
Ortho-ease handle
Adjust 30-39"
300lb weight capacity
\$35.50




Quad Cane, Std Grip
Adjust 30-39"
300lb weight capacity
6141A Small Base
6140A Large Base
\$35.50



Quad Cane, Bariatric
Adjust 29-38" 500lb Cap
Silver Vein finish
Contoured PVC grip
6130A Large Base
\$39.95

CANE, WALKER TIPS, ACCESSORIES

Cane Tips



Walker Tips



Inside Diameter/type	Color Available	Base Width	Price Each
1/2" Quad Cane Tip	Grey	1-1/4"	\$1.99
5/8" Standard Cane Tip	Grey or Black	1-1/8"	\$1.99
3/4" Standard Cane Tip	Grey or Black	1-1/4"	\$1.99
7/8" Standard Cane Tip	Grey or Tan	1-1/2"	\$1.99
5/8" Smooth Style Cane Tip	Black	1-1/2"	\$2.95
3/4" Smooth Style Cane Tip	Black	1-1/2"	\$2.95
1/2" Walker/Commode Tip	Black	7/8"	\$1.99
5/8" Walker/Commode Tip	Black or Grey	1"	\$1.99
3/4" Walker/Commode Tip	Black or Grey	15/16"	\$1.99
7/8" Walker/Commode Tip	Black or Grey	15/16"	\$1.99
1" Walker/Commode Tip	Grey	1-3/4"	\$1.99
1-1/8" Walker/Commode Tip	Grey	1-3/4"	\$2.20



Ice Gripper Attachment

5-point stainless steel, flips up when not in use 8320200, \$12.50



Elastic Cane Straps

One size fits all Brown or Black \$3.95



4-POINT TIP SELF-STANDING

CANE HOLDER

A70800

\$12.95 -EA



CANE CLIP UNIVERSAL

1188

\$9.95 -EA

ALUMINUM ADJUSTABLE CRUTCHES



Size	Item #	Fits Height Range	Price
Youth	EMP 4002	4'6" - 5'2"	\$35.00
Adult	EMP 4003	5'2" - 5'10"	\$35.00
X-Tall	EMP 4004	5'10" - 6'6"	\$35.00
XX-Tall	MIL MWD6500	5'7 - 6'10"	\$106.00

Supports Up To 300lbs

Lightweight, corrosion-resistant aluminum frame



Euro Style Forearm Crutches
MOV1014, \$25.99 -EA



Featherweight Lofstrand aluminum forearm crutch is made of aluminum casting for extra strength. Weight Capacity 250LBS \$130.00 -PR

Adult (Length)
Standard: 28-36"
Short: 24-32"
Tall: 33-41"

Child (Length)
Standard: 24-32"
Short: 20-27"



Quad Seat Cane, anodized aluminum, nylon seat webbing, weighs 2.4lbs, 250lb capacity
RTL10360, \$39.95