

Description	Add'l Info	Item #	Price
-------------	------------	--------	-------

DIAGNOSTIX 720 POCKET ANEROID SPHYG



- Precision crafted 300mmHg Japanese engineered manometer with classic black baked enamel finish exceeds industry performance standards
- Extended scaleplate graduation marks reduce parallax for easier viewing from virtually any angle
- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system prevents mis-cuffing
- Adflow™ filter screen check and deflation valve reduce dust build-up
- Adflow™ chrome-plated brass air release valve for precise deflation control
- Leatherette carrying case with durable nylon zipper
- Not made with natural rubber latex
- Assembled, inspected, and packaged in the U.S.A. from foreign components
- Meets or exceeds ANSI/AAMI/ISO 81060-1
- Lifetime warranty on gauge, 3 years on inflation system

\$65.50

720-11AAD..... ADIMALS	720-11ARB..... ROYAL BLUE
720-11ABCA..... BREAST CANCER AWARENESS	720-AGR..... GREEN
720-11ADG..... HUNTER GREEN	720-11AN..... NAVY
720-11AMT..... MEDICAL THEME	720-11AV..... PURPLE
720-11ABK..... BLACK	720-11AR..... RED

PROSPHYG 768 POCKET ANEROID SPHYG with MATCHING CASE



The 768 series features:

- Contemporary black enamel 300mmHg manometer provides day-in, day-out dependability exceeds industry standards
- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system, prevents mis-cuffing
- Chrome-plated brass air release valve
- Extra large, color-matched, dual pocket carry case, measures 10" x 6" x 3"
- Not made with natural rubber latex
- Assembled, inspected, and packaged in the U.S.A. from foreign components

768-11ARB..... ROYAL BLUE
768-11ADG..... HUNTER GREEN
768-11AN..... NAVY
768-11AV..... PURPLE
768-11ABK..... BLACK
768-11AR..... RED
768-11ABD..... BURGUNDY
768-11AOR..... ORANGE

\$55.95

AUTOMATIC BLOOD PRESSURE MONITORS



Product Features

- Clinically Validated for Accuracy
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- One Button Operation
- Store 60 Readings each for up to 4 users
- Average reading feature provides trending
- Time & date stamp for every reading
- AccuFit Plus Cuff fits most arms 8.6-16.5"
- Multi-user , upper armAC Power Adapter

PREMIUM WRIST BLOOD PRESSURE MONITOR

Advanced, multi-user blood pressure monitor that provides trusted accuracy and reliable measurements every time.

AND UB543, \$77.95

AND UA767FAC (WITH AC ADAPTER).....\$89.95





AND UA767F (WITHOUT AC ADAPTER).....\$79.95






DIGITAL BP MONITOR AUTO XL CUFF 16.5-23.6"	WITH MEMORY,AC ADAPTER	AND UA789AC	189.95
DIGITAL BP MONITOR AUTO SMALL ADULT CUFF	WITH MEMORY,AC ADAPTER	AND UA767PSAC	89.95



OTOSCOPE,POCKET XENON 3X MAG	INCLUDES:3-SPECS,2-AAA BAT	PRESTIGE 3470	65.00
PENLIGHT WITH PUPIL GAUGE DISP	EACH,6/PK	R89032	2.95

3470 XENON POCKET OTOSCOPE

Description	Add'l Info	Item #	Price
 <p style="text-align: center;">PULSE OXIMETERS</p> <ul style="list-style-type: none"> • 2 Year Warranty • Free Case Included! • Easy to Read Large LCD Display • FDA Approved for Accuracy • Efficient - Over 30 hours (1 AAA Battery Required) <p>DB12, \$89.95</p>	   <p>Fingertip measurement of SpO2% (Oxygen Saturation) and pulse rate. Measures quickly and accurately. LED display screen can display information in two directions. Low power consumption with auto shut-off feature. Up to 30 hours of continuous use. Lightweight and compact design.</p> <p>PRESTIGE 459-BLK BLACK PRESTIGE 459-HKW HOT PINK</p>		<p style="text-align: center;">\$39.95</p>





STETHOSCOPES											
 	<p style="text-align: center;">LITTMANN CARDIOLOGY - Diagnose and make decisions with confidence</p>										
	<ul style="list-style-type: none"> • Hear subtle changes in patient status • Useful in critical care and challenging environments such as the ED, ICU, Cardiac ICU, step-down unit and other dynamic locations • Ergonomic and high profile construction; 40% larger chestpiece and 60% deeper bell than the 3M™ Littmann® Classic III™ Stethoscope • Designed for use with adult and pediatric patients • Tunable, dual-sided stainless steel chestpiece with open or closed bell • Non-chill bell sleeve for greater patient comfort • Soft-sealing eartips provide optimal comfort and excellent sound occlusion • Made in the USA without natural rubber latex 										
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">STETHOSCOPE, LITTMANN CARDIOLOGY 4</td> <td style="width: 20%;">BLACK TUBE 27"</td> <td style="width: 15%;">3M 3128</td> <td style="width: 15%; text-align: right;">249.95</td> </tr> <tr> <td>STETHOSCOPE, LITTMANN MASTER CARDIOLOGY</td> <td>BLACK TUBE 27"</td> <td>3M 2160</td> <td style="text-align: right;">265.35</td> </tr> </table>				STETHOSCOPE, LITTMANN CARDIOLOGY 4	BLACK TUBE 27"	3M 3128	249.95	STETHOSCOPE, LITTMANN MASTER CARDIOLOGY	BLACK TUBE 27"	3M 2160	265.35
STETHOSCOPE, LITTMANN CARDIOLOGY 4	BLACK TUBE 27"	3M 3128	249.95								
STETHOSCOPE, LITTMANN MASTER CARDIOLOGY	BLACK TUBE 27"	3M 2160	265.35								
  	<p>Littmann's latest addition to their line. The Littmann Classic III features patented tunable technology that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece. High acoustic sensitivity for performing general physical assessments. The bell stays clear of dirt and debris by keeping it covered with the small diaphragm.</p>										
	STETHOSCOPE LITTMANN CLASSIC 3	BLACK TUBE 27"	3M 5620	126.50							
	STETHOSCOPE LITTMANN CLASSIC 3	NAVY TUBE 27"	3M 5622	126.50							
	STETHOSCOPE LITTMANN CLASSIC 3	CARIBBEAN BLUE TUBE 27"	3M 5623	126.50							
	STETHOSCOPE LITTMANN CLASSIC 3	BURGUNDY TUBE 27"	3M 5627	126.50							
	STETHOSCOPE LITTMANN CLASSIC 3	CIEL BLUE TUBE 27"	3M 5630	126.50							
	STETHOSCOPE LITTMANN CLASSIC 3	LAVENDER TUBE 27"	3M 5632	126.50							
	STETHOSCOPE LITTMANN LITEWEIGHT BLK	BLACK TUBE 28"	3M 2450	59.95							
	STETHOSCOPE LITTMANN LITEWEIGHT BRG	BURGUNDY TUBE 28"	3M 2451	59.95							
	STETHOSCOPE LITTMANN LITEWEIGHT CEIL	CEIL BLUE TUBE 28"	3M 2454	59.95							
	STETHOSCOPE LITTMANN LITEWEIGHT CRB	CARRIBEAN BLUE TUBE 28"	3M 2452	59.95							
	STETHOSCOPE LITTMANN LITEWEIGHT PINK	PINK TUBE 28"	3M 2456	59.95							
	STETHOSCOPE CLINICAL PLUS LTW BLK	31", 4.9oz SS STEM	PRESTIGE S123-BLK	39.95							
	STETHOSCOPE CLINICAL PLUS LTW PINK	31", 4.9oz SS STEM	PRESTIGE S123-HPK	39.95							
	STETHOSCOPE CLINICAL PLUS LTW NAVY	31", 4.9oz SS STEM	PRESTIGE S123-NAV	39.95							
STETHOSCOPE CLINICAL PLUS LTW PURPLE	31", 4.9oz SS STEM	PRESTIGE S123-PUR	39.95								

The Adtemp™ 418N features:

- 8-second measurement using proprietary predictive technology
- Range 93.2°F to 111.0°F ±.2°F or 34.0°C to 43.9°C ±.1°C depending upon scale selection
- Large LCD display with color-coded backlighting (green for readings in the normal range and red for readings above 99.5°F/37.5°C)
- Replaceable 3V battery (CR2032) provides up to 2,000 measurements
- Measures 5.31"L x 1.14"W x 0.59"H, weighs 0.77 oz (22 g) with battery
- Integral scabbard-type carry case, Includes five probe sheaths
- Meets the accuracy requirements for intermittent thermometers under ASTM E1112

ADC 418N, \$13.99



Description	Add'l Info	Item #	Price
NON-CONTACT THERMOMETERS			
	<p>The Adtemp™ 429 Non-Contact Thermometer features an infrared sensor that reads temperature from the middle of the forehead with no skin contact required for a safe and comfortable experience that minimizes the risk of cross-contamination.</p> <p>The Adtemp™ 429 Infrared Thermometer features:</p> <ul style="list-style-type: none"> • Infrared sensor simplifies temperature measurement—ideal for screening • Measures in about 1 second without patient contact eliminating the concern of cross contamination or costly probe covers • Measurement feedback system signals proper operating distance for accuracy • Measures at a distance of 4-6cm (1.5"-2.5") from patient's forehead • Dual mode - body and object • Large backlit LCD display is easy to read • Wide Temp Range • 30 reading memory • Fahrenheit-Celsius switchable • Auto-off function after 30 seconds conserves battery life • Audible beep when measurement is complete, Night mode mutes volume • Compact ergonomic design 5" x 2.06" x 1.75" (129.9mm x 52.4mm x 44.5mm), weighs approximately 3oz (87g) with batteries • Optional cradle available for wall mounting <p>Uses 2 AAA batteries (included)</p> <p>ADC 429, \$ 79.99 2 year warranty</p>		
	<p>Adtemp Mini 432 Non-Contact Thermometer</p> <p>ADC's newest non-contact thermometer, the Adtemp Mini 432, is compact, lightweight, and fast. Offering an ingenious design and incredible portability – along with the safety and convenience of a non-contact instrument – the 432 is sure to become an essential instrument for busy professionals.</p> <ul style="list-style-type: none"> • One-second readings • Compact, lightweight size, measures 3.5" and weighs 1.6 oz • No skin contact required -- minimize the risk of cross contamination • Color-coded LCD display with fever warning • 9-reading memory, °F/°C switchable display • Rechargeable Li-ion battery, micro USB cable included • 1,000 measurements on a single charge, 300,000 measurements lifetime • Clips to wearable, retractable lanyard (included) so it's always on hand • Audible beep when measurement is complete <p>ADC 432, \$69.95 2-year warranty</p>		
 <ul style="list-style-type: none"> • Non-Contact, Fast , Accurate Readings • Measures Body, Surface and Ambient Temperatures • Backlit Digital LCD Display Screen • High Temperature Alarm • Measures in Both Fahrenheit & Celsius • Features 32 data memory <p>Includes two AAA batteries, 1 year warranty</p> <p>DT-29, \$59.95</p>		 <ul style="list-style-type: none"> • 2 modes, 1 push: Body & object mode • 1-second accurate measurement time • Auto turn-off function for energy saving • 32 readings memory recall 1 push switch • Fahrenheit °F or Celsius °C <p>SM017011CF, \$49.95 -Each</p>	