

 <p>VERSATILT OVERBED TABLE 15x30" WALNUT MELANINE TOP 26.5-39" INFINITE HT ADJUST 950H, \$217.00</p>	 <p>FOOT STOOL, 300LB CAP 10550, \$32.50</p>	 <p>FOOT STOOL WITH HANDLE 10552, \$55.50</p>	 <p>FOOT STOOL HD 500LB CAPACITY 10551, \$45.50</p>
---	--	--	---

PATIENT LIFTS



- Silver Vein steel construction
 - 6-point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
 - High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
 - Choose either 3" or 5" casters
 - Adjustable-width base
 - 450 LB weight capacity
 - Lifetime warranty (frame), 1 year limited (pump)
- 13023SV, \$695.00

The PL400HE is Bestcare's answer to home care needs. Featuring a 400 LB maximum weight capacity and Bestcare battery powered control system for power lifting, the PL400HE is easily adaptable for a wide range of challenges while still maintaining a high standard of comfort for patients. The use of Bestcare's battery powered electronic system provides a low-cost solution for power lifting while providing convenience for caregivers. Designed for home care use and effective at providing quality care.



- 400 LB carrying capacity
 - Six-point spreader bar
 - Electronic control system
 - Optional digital weigh scale
 - Optional two-point spreader bar
- PL400HE, \$1,895.00



MESH SLING FULL BODY (NO COMMODE OPEN), 450LB CAPACITY

SL-R110 MEDIUM: FITS WEIGHT RANGE 100-195 LBS, \$95.00
SL-R111 LARGE: FITS WEIGHT RANGE 150-300 LBS, \$96.00



MESH SLING FULL BODY (COMMODOE OPEN), 450LB CAPACITY

SL-R114 MEDIUM: FITS WEIGHT RANGE 100-195 LBS, \$96.00
SL-R115 LARGE: FITS WEIGHT RANGE 150-300 LBS, \$98.00



VERSA-HELPER TRAPEZE

- Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance
 - Rubber bumper guards protect walls
 - Hand bar location fully adjustable
 - Height adjustable from 43" to 52"
 - VersaGuard coated clamps and hand bar
 - For use with metal frame headboards only
 - Assisted weight capacity of 250lbs
- 2800GA (TRAPEZE ONLY), \$199.95
2840GA, (FLOOR BASE), \$125.00



Floor Base

ELECTRIC HOSPITAL BEDS



- Constructed from durable steel with a spring deck and split-frame design
- Spring deck is adjustable from 15.7"-22.8"
- Motor is compatible with a 9-volt battery to adjust bed position in a power outage
- Weight capacity 450 LBS
- Semi-Electric beds: adjustable head and foot sections, overall height it manually operated with a crank
- Full Electric Beds: all 3 functions (head, foot, overall height) are operated with a push-button pendant

HOSPITAL BED RAILS



HALF RAILS

- Unique design prevents individuals from being lodged between bars
- Ideal for maximum protection when bed is elevated
- Fits most hospital beds with steel frame
- Attractive brown-vein finish is easy-to-maintain and less institutional looking
- Constructed of 1" steel



FULL RAILS

- Attractive brown-vein finish
- Constructed of 1" steel
- Telescopes with movement of bed
- Easily adjusts up or down with spring-loaded release
- Spring-loaded crossbars are easy-to-install and remove

ACCESSORY ITEMS

HOSPITAL BED FITTED SHEET 36X80'

50-50 COTTON POLYESTER \$22.25 -EA

THERMAL BLANKET 100% COTTON

66x90" \$25.50 -EA

HOSPITAL BED FITTED SHEET SET

1-FITTED,1-FLAT,1-PILLOWCASE
50-50 COTTON POLYESTER \$45.95 -EA



BLANKET SUPPORT

FITS ALL BEDS, ADJUSTS TO 23" HIGH
ALUM/BRONZE FINISH
KEEPS BLANKET OFF FEET \$29.95 -EA

BED RAIL BUMPER PADS

72Lx14Hx1.5" \$69.95 -PR


THERAPEUTIC 5-ZONE SUPPORT MATTRESS
 Deluxe horizontal, cross-cut foam mattress provides comfort, support, and pressure redistribution over 5 therapeutic pressure zones



- Better compliance and durability than innerspring mattresses allow for improved patient satisfaction and fewer service calls
 - Nonskid bottom keeps mattress safely in place
 - Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633
 15019, \$329.95

GRAVITY-7 MULTI-LAYERED MULTI-ZONED MATTRESS

High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction



- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
 - Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
 - Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
 - High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
 - Visco foam heel section offloads pressure from vulnerable heels
 - Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633
 15770, \$379.95


Air-Pro™ Series

Alternating Pressure Pad & Pump Overlay Mattress Systems

AIR-PRO™ SERIES - OVERLAY MATTRESS SYSTEM

The Air-Pro™ Plus & Elite Alternating Pump and Overlay Mattress Systems include a foldaway bed hanger bracket. This feature keeps the pump off the floor and reduces any accidental damage. Our easy access filter located on the bottom of the pump reduces any chance of dust particles from entering the pump compressor motor. We include a heavy duty, deluxe bubble style pad with head and foot tuck under flaps. The adjustability of the Elite allows this system to be used for heavier patients. Air-Pro™ meets California Technical Bulletin #117 for fire retardancy.

Shown with 4810 pad



AIR-PRO SERIES ALTERNATING PRESSURE OVERLAY SYSTEM

4400SCB Complete System, \$129.95
 4800 Pad only, \$29.95



Model 4800 Pad

AIR-PRO™ SERIES - OVERLAY MATTRESS PAD

Material: Medical grade PVC with Anti-microbial, fire-retardant treatments.

Design: Bubble Style with more touch points for effective pressure redistribution.

MED-AIRE 8" ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM

Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control.
 36Wx80Lx8H" Fluid resistant nylon cover



- Designed to prevent, treat and heal pressure ulcers in the home or long-term care setting
- Control unit offers 10-minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Audio and Visual alarm alerts when air pressure is low

- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
 - 8" deep air cells are easily removed and replaced with air leakage plugs
 - Removable Nylon stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
 - CPR valve allows for rapid deflation
 - Tubing quickly and easily disconnects with EZ Lock/Quick Release design
- 14027, \$1299.00

Premium Guard Gel Overlay

Best in class gel overlay optimizes pressure redistribution and is now available in up to 500lbs weight capacity

- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length "head-to-toe" gel bladders minimize "bottoming out" and "pitting"
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration



14893 PremGuardGel/FoamOverlay34x76x3.5" 300lb Capacity, \$209.95

14901 PremGuardGel/FoamOverlay42x76x3.5" 500lb Capacity, \$495.00

HOSPITAL BED PACKAGE PRICING, INCLUDES DELIVERY AND SET-UP*				
PKG#	Bed Type	Rail Type	Mattress	Price
SF5Z	SEMI-ELECTRIC	FULL RAIL	15019 MATTRESS	\$1,299.00
SFG7	SEMI-ELECTRIC	FULL RAIL	15770 MATTRESS	\$1,350.00
SH5Z	SEMI-ELECTRIC	HALF RAIL	15019 MATTRESS	\$1,299.00
SHG7	SEMI-ELECTRIC	HALF RAIL	15770 MATTRESS	\$1,350.00
FF5Z	FULL-ELECTRIC	FULL RAIL	15019 MATTRESS	\$1,499.00
FFG7	FULL-ELECTRIC	FULL RAIL	15770 MATTRESS	\$1,550.00
FH5Z	FULL-ELECTRIC	HALF RAIL	15019 MATTRESS	\$1,499.00
FHG7	FULL-ELECTRIC	HALF RAIL	15770 MATTRESS	\$1,550.00
*within 20-mile service area, delivery surcharge added if over				

